

2915 Everest Lane  
 Plymouth, MN 55447  
 www.eyesurgin.com

(763) 473-2533  
 Fax:(763) 208-8105



**Prices quoted FOB Plymouth, MN, exclusive shipping. Terms: Net 30. V/MC accepted.**

<b>CLEARSCAN® UBM and B Scan Sterile Probe Cover</b> 1 Box (10 units)	\$95.00
<b>Prager Shell® Package I:</b> 1 Prager Shell, 25 Tubing Kit B, 1 Cleaning/Storage Tray	\$285.00
<b>Prager Shell® Package II:</b> 1 Prager Shell, 25 Tubing Kit B	\$265.00

**ESI PRAGER SHELL MODELS:**

<u>Model Number</u>	<u>Ultrasound Companies</u>
P-100	Accutome / Micro Medical Devices
P-100-OCU	Ocuserv
P-101	DGH 2009 and earlier model
P-101-2	DGH 2010 and later model
P-102	Tomey
P-103	Sonomed / Ellex / OPKO
P-103 SV	Sonomed / Ellex / OPKO 12.5mm diameter
P-104	Quantel Medical
P-104 SV	Quantel Medical 12.5mm diameter
P-107-2	Nidek
P-105	Alcon OcuScan
P-110 RXP	Alcon RXP OcuScan and UltraScan A/B
P-110 RXP SV	Alcon RXP OcuScan and UltraScan A/B 12.5mm diameter

**Prager Shell Autoclavable Gold Series available by special order.**

**Kit B**

Designed for infusing BSS from a squeeze bottle into the Prager Shell®. Each end of the tubing has a different type of Luer connector. **Luer connector features a one-way flow check valve to prevent back flow of** contaminated fluid into the BSS source. Kit B is for single patient use.

**Kit B** Product number: B1263, 1 Box (25 units per box) \$58.75

**Kit A**

**Does not have the one-way flow check valve** for preventing back flow of patient contaminated fluid into the BSS source. Kit A should only be used with a syringe. Kit A is for single patient use.

**Kit A** Product number: B1236, 1 Box (25 units per box) \$52.50

**Storage/Cleaning Tray:**

Custom tray holds shell and probe for conveniently soaking in an aseptic solution after patient use. Protects the Prager Shell® while stored.

Tray Product number: T200

1 Tray \$44.00

1 Tray with an initial order of 50 or more Kit A/Kit B. One time offer per location. \$22.00