

ESI, Inc.

Eye Surgical Instruments

2915 Everest Lane North
Plymouth, MN 55447

(763) 473-2533
(763) 208-8105 Fax
sales@eyesurgin.com



Prices listed are effective 2/1/18. FOB Plymouth, MN, Terms: Net 30.

Check or V/MC/AMEX accepted.

ClearScan® CS200: UBM and B Scan Sterile Probe Cover 1 Box (10 units)	\$105.00
Prager Shell® Package I: 1 Prager Shell®, 25 Tubing Kit B, 1 Cleaning/Storage Tray	\$315.00
Prager Shell® Package II: 1 Prager Shell®, 25 Tubing Kit B	\$295.00

Prager Shell® Models:

Model Number	Product Description
P-100	Accutome / MicroMedical
P-100AS	Accutome Probe Auto Stop 14.9mm diameter
P-101	DGH 2009 and earlier model 14.9mm diameter
P-101-2	DGH 2010 and later model 14.9mm diameter
P-102	Tomey AL100 14.9mm diameter
P-102-2	Tomey AL 4000 14.9mm diameter
P-103-2	Sonomed 14.9mm diameter
P-103-2GS	Sonomed 14.9mm Autoclavable Gold Series
P-103-2SV	Sonomed 12.5mm diameter
P-104	Quantel Medical 14.9mm diameter
P-104GS	Quantel Medical 14.9mm Autoclavable Gold Series
P-104SV	Quantel Medical 12.5mm diameter
P-105	Alcon Ocuscan 14.9mm diameter
P-106	Ellex 14.9mm diameter
P-106GS	Ellex 14.9mm Autoclavable Gold Series
P-106SV	Ellex 12.5mm diameter
P-107-2	Nidek 14mm diameter
P-110RXP	Alcon RXP OcuScan and UltraScan A/B 14.9mm diameter
P-110RXPVS	Alcon RXP OcuScan and UltraScan A/B 12.5mm diameter

Prager Shell® Autoclavable Gold Series available by special order and \$40.00 surcharge.

Kit B: Product Number 011-B1263

Designed for infusing BSS from a squeeze bottle into the Prager Shell®. **Luer connector features a one-way flow check valve** to prevent back flow of patient contaminated fluid into the BSS source. Kit B is for single patient use.

Kit B 1 Box (25 units per box) \$62.50

Kit A: Product Number 011-B1236

The one-way flow check valve is not included for preventing back flow of patient contaminated fluid into the BSS source. Kit A should only be used with a syringe. Kit A is for single patient use.

Kit A 1 Box (25 units per box) \$58.75

Storage/Cleaning Tray: Product Number T-200

Custom tray holds shell and probe for conveniently soaking in an aseptic solution after patient use. Protects the Prager Shell® while stored.

T200 Tray \$44.00